TITRATION VALIDATION TEST

NAME:			25 Marks + 5 for accuracy = 30 marks		
				Mark =	/30
AIM To determine the p	ercentage by mass of am	monia in house	ehold ammonia		
METHOD	ana ant with the annuarie		معدد المعالمة	:f	
2. Pipette 20 mL o	oment with the appropria of household ammonia in monia solution. Do not pi	ito a 250 mL vo	lumetric flask. I	Determine the ex	
Make the volur doing this).	ne up to 250 mL with dist	tilled water and	d mix well. (Use	e the appropriate	technique for
4. Titrate the star	ndard HCl solution from a with 3-5 drops of methyl o sults below.	_		s of diluted house	ehold ammonia
6. Repeat steps 3	and 4 until you obtain 3 o	concordant res	sults.		
RESULTS					
Concentration of the standard HCl solution =					
Mass of the 250 mL volumetric flask =					
Mass of the 250 mL volumetric flask + 20mL of household ammonia =					
Mass of 20 mL of household ammonia =					

	Rough estimate	Titration 1	Titration 2	Titration 3	Average Titration volume (mL)
Initial reading (mL)					
Final reading (mL)					
Titration volume (mL)					

(1 mark)

ANALYSIS OF RESULTS

 Write an appropriate equation for the reaction between the ammonia and hydrochloric a calculate the percentage by mass of ammonia in the household ammonia. 		
	-	
_		
	_	

CCGS Year 12 Chemistry			Acid -base Titration Validation Test 2019		
				(7 marks)	
Accura	cy: Mark =/5				
QUESTI	ONS				
1.(a) C	omplete the following table b	y circling the corre	t alternative in the right hand	column	
	Equipment		Rinsed with		
	Volumetric flask	Water	Ammonia solution		
	Burette	Water	Hydrochloric acid		
	Conical flask	Water	Ammonia solution		
• •	ow would the titration volumen part (a)? Justify your answer.		burette was rinsed with the in	(3 marks) correct solution	
				(2 marks)	
ri w	nse it with ammonia solution l	before titrating. As	I their pipette with water but r suming they do only one titrat mmonia in their household am	ion, what effect	

(4 marks)

(b) Use your graph to explain why methyl orange (pH range 3.2 – 4.4) was chosen as the indicator in preference to phenolphthalein (pH range 8.2 – 10) for this titration. Use equations to support your answer.

CCGS Year 12 Chemistry	Acid -base Titration Validation Test 2019		
			

END OF TEST